



*Variable Volume Fully Autoclavable
Pipettes VVP109L*

www.labodam.com // info@labdam.com

Description :

Specifications :

Volume Range	1000 - 5000 μL
Increment	50 μL
Test Volume	5000 μL
	2500 μL
	1000 μL
Accuracy error	0.50 %
	0.60 %
	0.70 %
Precision error	0.15 %
	0.30 %
	0.30 %

Features :

Autoclavable at 121 $^{\circ}\text{C}$

Material: Polycarbonate

Available in 0.1 μL to 5000 μL volume

Ergonomic design provides excellent operating experience

Easy-to-read volume display

Trouble-free calibration and maintenance



Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA

United Kingdom

www.labodam.com // info@labodam.com