



*Variable Volume Multi Channel Fully  
Autoclavable Pipette VMP100L*

## Description :

## Specifications :

No of Channels	8
Volume Range	0.5 -10 $\mu\text{L}$
Increment	0.1 $\mu\text{L}$
Test Volume	10 $\mu\text{L}$
	5 $\mu\text{L}$
	1 $\mu\text{L}$
Accuracy error	1.50 %
	2.50 %
	4.00 %
Precision error	1.50 %
	2.50 %
	4.00 %

## Features :

Autoclavable at 121  $^{\circ}\text{C}$

Material: Polycarbonate

Volume settings can be clearly read through digital display

Effortless pipetting accessibility with rotatable Dispensing head

Easy repair and maintenance with separate piston and tip cone assembly

Compatible with most universal tip



**Labodam Equipment Ltd**

**18a Melton Road Leicester LE4 5EA**

**United Kingdom**

**[www.labodam.com](http://www.labodam.com) // [info@labodam.com](mailto:info@labodam.com)**