



*Variable Volume Multi Channel Fully
Autoclavable Pipette VMP101L*

Description :

Specifications :

No of Channels	8
Volume Range	5-50 µL
Increment	0.5 µL
Test Volume	50 µL
	25 µL
	5 µL
Accuracy error	1.00 %
	1.50 %
	3.00 %
Precision error	0.50 %
	1.00 %
	2.00 %

Features :

Autoclavable at 121 °C

Material: Polycarbonate

Volume settings can be clearly read through digital display

Effortless pipetting accessibility with rotatable Dispensing head

Easy repair and maintenance with separate piston and tip cone assembly

Compatible with most universal tip



Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA

United Kingdom

www.labodam.com // info@labodam.com