



*Variable Volume Multi Channel Fully
Autoclavable Pipette VMP102L*

Description :

Specifications :

No of Channels	8
Volume Range	50-300 µL
Increment	5 µL
Test Volume	300 µL
	150 µL
	50 µL
Accuracy error	0.70 %
	1.00 %
	1.50 %
Precision error	0.25 %
	0.50 %
	0.80 %

Features :

Autoclavable at 121 Â°C

Material: Polycarbonate

Volume settings can be clearly read through digital display

Effortless pipetting accessibility with rotatable Dispensing head

Easy repair and maintenance with separate piston and tip cone assembly

Compatible with most universal tip



Labodam Equipment Ltd

**18a Melton Road Leicester LE4 5EA
United Kingdom**

www.labodam.com // info@labodam.com