



*Variable Volume Multi Channel Fully  
Autoclavable Pipette VMP105L*

## Description :

## Specifications :

No of Channels	12
Volume Range	50-300 $\mu\text{L}$
Increment	5 $\mu\text{L}$
Test Volume	300 $\mu\text{L}$
	150 $\mu\text{L}$
	50 $\mu\text{L}$
Accuracy error	0.70 %
	1.00 %
	1.50 %
Precision error	0.25 %
	0.50 %
	0.80 %

## Features :

Autoclavable at 121  $\hat{\text{A}}^{\circ}\text{C}$

Material: Polycarbonate

Volume settings can be clearly read through digital display

Effortless pipetting accessibility with rotatable Dispensing head

Easy repair and maintenance with separate piston and tip cone assembly

Compatible with most universal tip



**Labodam Equipment Ltd**

18a Melton Road Leicester LE4 5EA

United Kingdom

[www.labodam.com](http://www.labodam.com) // [info@labodam.com](mailto:info@labodam.com)