



*Variable Volume Multi Channel
Pipette VMP200L*

Description :

Specifications :

No of Channels	8
Volume Range	0.5 -10 μL
Increment	0.1 μL
Test Volume	10 μL
	5 μL
	1 μL
Accuracy error	1.50 %
	2.50 %
	4.00 %
Precision error	1.50 %
	2.50 %
	4.00 %

Features :

Autoclavable at 121 $^{\circ}\text{C}$

Material: Polycarbonate

Volume settings can be clearly read through digital display

Effortless pipetting accessibility with rotatable Dispensing head

Easy repair and maintenance with separate piston and tip cone assembly

Compatible with most universal tip

labodam

Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA

United Kingdom

www.labodam.com // info@labodam.com