



Variable Volume Pipette VVP210L

Description :

Specifications :

Volume Range	2 - 10 mL
Increment	0.1 mL
Test Volume	10 mL
	5 mL
	2 mL
Accuracy error	0.60 %
	1.20 %
	3.00 %
Precision error	0.20 %
	0.30 %
	0.60 %

Features :

Autoclavable at 121 °C

Material: Polycarbonate

Available in 0.1 µL to 10000 µL volume

Ergonomic design provides excellent operating experience

Easy-to-read volume display

Trouble-free calibration and maintenance



Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA
United Kingdom

www.labodam.com // info@labodam.com