



Rapid Low Temperature Circulating Bath LD-LRTB-A12

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Rapid Low Temperature Circulating Bath LD-LRTB-A12

Overview

Rapid Low Temperature Circulating Bath LRTB-A12 is a microcomputer controlled 3 L bath with cold ring liquid circulation and cryogenic liquid with a temperature range of -10°C . Refrigeration unit with low temperature cooling pump facilitates consistent high pressure and maintenance of low temperature under working conditions of reaction experiments. Equipped with anti-corrosive circulatory system preventing traffic flow, high pressure and low pump body temperature.

Features:

- Volume - 3 L
- Low temperature cooling pump
- Microcomputer controlled temperature
- Temperature and time LCD display
- Built-in liquid level monitor
- Refrigeration unit with relay , protector and capacitor
- Anti-corrosive circulatory system

Specifications:

Volume	3 L
Pump flow	20 L / min
Head of delivery	10 m
Temperature range	- 20 °C
Temperature stability	$\pm 0.5^{\circ}\text{C}$
Refrigerating capacity	782 W