



Shaking Incubator (Stack Type) LD-LSI-D10

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Shaking Incubator (Stack Type) LD-LSI-D10

Overview

Shaking Incubator (Stack Type) LSI-D10 has a temperature range from 4 °C to 60 °C and its rotating speed is from 50 to 300 rpm. The vibrating range of the device is ± 25 mm. It stand with a maximum capacity of 2000 ml — 12 or 1000 ml — 18 or 500 ml — 35 or 250 ml — 40 with a temperature accuracy of ± 0.1 °C. The device uses microprocessor PID control system for smooth functioning. Additionally it has an automatic alarm system for over-temperature, over-speed and for over-current to protect the device. It also consist of advanced Unishaft drive with low noise operation. The inner wall is made of stainless steel and has good corrosion resistance. And the external wall is made of cold-rolled steel with anti-bacteria power coating.

Features:

- Advanced Unishaft drive with low noise operation
- Alarm function is available for audio as well as visual
- When the motor is overheating and temperature is out of control automatic power-off protection system is activated
- Shaking plate of the device is made of electroplating aluminum, anti-rust or anti-corrosion
- Eight self-compiled programs, with different speed and time setting
- Automatic operation, auto-stop, timing, time display, parameters memory and recovery function
- LCD display is used to show the result of the experiment

Applications:

Shaking Incubator (stack type) is used for the studies of solubility. It can be also used for growth of any kind of cell including bacterial cultures, tissue cultures and yeast incubator.

Specifications:

Display	LCD
Drive Mode	Unishaft drive
Consumption	1120 W
Refrigerant	CFC-Free, 134 A
Gross Weight	410 kg (per unit)
Package Size	1400×1010×1130 mm (per unit)
Power Supply	AC 110 V / 220 V $\pm$ 10 %, 50 / 60 Hz
Shaking Mode	Rotatory Shaking
Timing Range	0 to 500 hours
External Size	1300×910×980 mm (per unit)

Internal Size	870×550×395 (per unit) mm
Rotating Speed	50 to 300 rpm
Safety function	Audio and visual alarm for over-temperature; Alarm for over-speeds; Independent over-temperature protector; Shutdown protection for door opening; Leakage or over-current protection
Vibration Range	± 25 mm
Circulation Mode	Forced convection
Maximum Capacity	2000 ml × 12 or 1000 ml × 18 or 500 ml × 35 or 250 ml × 40
External Material	Cold-rolled steel with anti-bacteria power coating
Internal Material	Mirror stainless steel
Rotating Accuracy	± 1 rpm
Standard Capacity	250ml × 40
Temperature Range	4 °C to 60 °C
Shaking Plate Size	850×560 mm
Cooling Consumption	200 W
Heating Consumption	800 W
Temperature Accuracy	± 0.1 °C
Temperature Uniformity	± 0.25 °C