



Crude Petroleum Water Content Tester LD-LWT-B10

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Crude Petroleum Water Content Tester LD-LWT-B10

Overview

LWT-B10 is crafted to determine water content in crude petroleum according to ASTM D4006 “Test Methods for Water Content of Crude Petroleum” international standard.

Features:

- Simple Design
- Adjustable Electric Heating Power
- Adopts Heating Mantle

Applications:

It is used to determine water content in crude oil, refined petroleum products, fuels, biofuels, lubricants, chemicals, and other products

Specifications:

Condenser	400 mm ± 5 mm
Dimensions	820 × 500 × 260 mm
Power Supply	AC 220 V ± 10 %, 50 Hz
Water receiver	5 ml
scale division	0.05 ml
Relative Humidity	≥ 85 %
Distillation flask	1000 ml
Ambient Temperature	(-10) °C ~ 35 °C