



## **Corrosive Cabinet LD-LCOC-B12**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Corrosive Cabinet LD-LCOC-B12

## Overview

Corrosive Cabinet LCOC-B12 is a dual walled manual combustion cabinet with Capacity 83 L and 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features three adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

## Features:

- Capacity - 83 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Three adjustable shelves
- Adjustable spill retaining shelves

## Applications:

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

## Specifications:

|                          |                     |
|--------------------------|---------------------|
| Weight                   | 101 kgs             |
| Capacity                 | 83 L                |
| Door type                | Single door manual  |
| Adjustable shelf         | 3                   |
| Dimension ( H x W x D )  | 1650 x 600 x 460 mm |
| Shelves loading capacity | 50 kgs              |