



Paraffin Block Cabinet LD-LPBC-A10

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Paraffin Block Cabinet LD-LPBC-A10

Overview

Paraffin Block Cabinet is metallic standalone unit, comes with 18 drawers. Paraffin cabinet is used to store paraffin blocks, Drawers is consisting of ABS slide plate and cabinet material is cold rolled steel plate and surface is treated with Degree oil, Table regulator, Zinc Parkerizing, Passivation, Powder spraying which is resistant to most chemical aggression.

Features:

- 18 Drawers with loading capacity of 12000 paraffin block
- Slot is made by special mold
- Corrosion resistant material

Applications:

Paraffin block cabinet is used for efficient and organized storage of paraffin blocks across histology, pathology, forensics, research centers etc.

Specifications:

Capacity	18 Drawers (3 layers x 6 grid, plus a base) Loading: 12000 pcs paraffin block
Mark groove	Integration forming
Gross weight	29kg, 3pcs; 33kg, 1pc
Slot material	Made by special mold
Drawer material	Specialized slide drawer, ABS slide plate
Cabinet material	Cold-rolled steel plate. Surface Treatment: Degree oil, Table regulator, Zinc Parkerizing, Passivation, Powder spraying
Package Size (W x D x H)	530 x 570 x 440 mm, 3pcs; 530 x 570 x 565 mm, 1pc
External Size (W x D x H)	450 x 478 x 1295 mm