



Standard strip cutter LD-TP-G10

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Standard strip cutter LD-TP-G10

Overview

Standard strip cutter TP-G10 with sleek and advanced design is used to cut the sheet of cardboard or paper into fine strips. It's a common standard cutter which can cut the sample with 15 mm width. With user-friendly operation, it finds applications in paper & printing industry.

Features:

- Advanced and sleek design
- Sample width 15 mm
- Sample maximum length 300 mm
- Easy to use

Applications:

Strip cutter is used in paper, packaging, inspection department, for paper crafts making industry.

Specifications:

Weight	4 kg
Dimension	240×520×70 mm
Parallelism	± 0.1 mm
Sample width	15 mm
Sample maximum length	300 mm