




## **Polypropylene Cell Spreader LD-PCS300L**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Polypropylene Cell Spreader LD-PCS300L

## Features:

- Material: Polypropylene
- Ideal for bacterial CFU counts
- Gamma ray sterilized
- DNase, RNase and pyrogenic free

## Specifications:

|                   |                      |
|-------------------|----------------------|
| Size              | 14.5 x 3.8 mm        |
| Color             | Yellow               |
| Shape             | L                    |
| Packaging         | 1 pcs/bag            |
| Sterilization     | Gamma ray sterilized |
| Quantity per case | 1000                 |