



## **Racked Bar-Coded Cryogenic Vials LD-BCV107L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Racked Bar-Coded Cryogenic Vials LD-BCV107L

## Features:

- Material: Polypropylene
- Low Temperature Storage: -196 Â°C to 121 Â°C
- Sterilization Level: SAL 10-6
- Chemically resistant and scratch-proof 2D Laser Etched code on Bottom
- U-shaped for maximum working volume
- Barcoded rack for tracking
- Inner Diameter: 0.7 mm
- Cap Color: Transparent, White, Blue, Orange and Green
- DNase, RNase and Pyrogenic free

## Specifications:

|                   |                            |
|-------------------|----------------------------|
| Height            | 28.1 mm                    |
| Barcode           | 1D                         |
| Cap Type          | Push Cap                   |
| Capacity          | 1.4 mL                     |
| Packaging         | 96 pcs/Rack, 48 Racks/Case |
| Working Volume    | 1.19 mL                    |
| Quantity per Case | 4608                       |