



## **Ultra Low retention Pipette Tips LD- LRT101L**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Ultra Low retention Pipette Tips LD-LRT101L

## Features:

- Universal compatibility
- Material: Polypropylene
- Autoclavable at +121°C
- DNase, RNase and Human DNA free
- Micro Point Design
- Pre Sterilized
- Long reach design
- Optional Color “ Yellow (LRT106L, LRT107L, LRT108L, LRT109L, LRT110L and LRT111L)
- Optional Color “ Blue (LRT112L, LRT113L, LRT114L, LRT115L, LRT116L and LRT117L)

## Specifications:

|                   |                             |
|-------------------|-----------------------------|
| Type              | Long                        |
| Capacity          | 10 µL                       |
| Packaging         | 1000 Tips/Bag, 10 Bags/Case |
| Length of tube    | 44.60 mm                    |
| Quantity per Case | 10000                       |